

Metkon Application Note

SAMPLE: Preparation of Hardened Gear sample

INDEX

1. INTRODUCTION

2. APPLICATION REQUIREMENTS

**3. SAMPLE PREPARATION
PROCESSES**

4. RESULT

In this application, hardened gear sample used which is shown in the picture below.

The purpose of the investigation is to describe hardening layer and microstructure.



Sample dimensions and requested cutting directions. (Hardness 45~50 HRC)

APPLICATION REQUIREMENTS

A. SECTIONING



	Order Code	Description
Equipment Used	14 56	Servocut 301 - MA
Clamping Device	15 01 GR 0170	MBU 1031 Quick clamping Device with clamping shoe Vise assembly, Left
Cutting Fluid	19-902	Metcool, Nature Friendly Soluble Oil, 5lt.
Cutting Disc	19-042	TRENO-M Ø 300 mm, for Medium Hard Steels > 23-55 HRC <

B. MOUNTING




	Order Code	Description
Equipment Used	25 07	Ecompress 100
Mould Assembly	26 06 - 02	Mould assembly 40 mm diameter
Mounting Powder	29-012 29-011	Diallyphtalat, 1 kg. Epoxy, Hard 1 kg.

C. GRINDING & POLISHING





	Order Code	Description
Equipment Used	45 03	DIGIPREP 251 Grinding & Polishing System
Equipment	31 22	Aluminum wheel, 250 mm
Accessories	31 63	Splash Guard, 250 mm
	39-003-250	Ø 250 mm, Special Magnetic Foil
	39-093-250	Ø 250 mm, Thin Metal Plate(5 pcs.)
Sample Holder	33 01	Specimen holder, 6 x Ø40 mm

SAMPLE PREPARATION PROCESSES



The sample is attached as it shown in the above photo with the **15 01** clamping device.



After first cutting the sample moved a little to get a slice and get gear teeth

The cutting parameters are below;



Table feed rate is adjusted to; **200 μ / sec**

Rpm is adjusted to; **2200**

Travel is adjusted to; **55 mm**

Force is adjusted to; **7A**

With these parameters, 2 cutting process took approx. 6-7 minutes.



The final piece divided into two parts for molding operation





ECOPRESS 100 was used for mounting operation with below parameters;

Heating temperature: 180 °C

Pressure: 240 bar

Heating time: 3 min

Cooling type: Standard cooling

Cooling temperature: 35 °C

The Mounting process took approx. 10 minutes,


We used Digiprep 251 machine to grind and polish the samples with following parameters.

The grinding & polishing process took app. 10 min.

	<i>Surface</i>	<i>Abrasive</i>	<i>Lubricant</i>	<i>Force per sample(N)</i>	<i>Time(min.)</i>	<i>Disk speed(rpm)</i>	<i>Head speed(rpm)</i>
Grind. Step 1	MAGNETO I [38-040-054]	54 micron Diamond	Water	20 N	1 min.	200	100
Grind. Step 2	MAGNETO II [38-040-018]	18 micron Diamond	Water	25 N	2 min.	250	100
Polishing Step 1	SATO cloth [39-013-250]	6 μ diamond susp. [39-430-M]	DIAPAT [39-502]	25N	3 min.	200	100
Polishing Step 2	FEDO-1 cloth [39-065-250]	1 μ diamond susp. [39-410-M]	DIAPAT [39-502]	20N	3 min.	200	100
Etching	3% Nital solution						



After grinding & polishing operation



After etching procedure (hardened layer can be seen)

As a result hardened sample was subjected to the following operations;

Cutting → Mounting → Grinding → Polishing → Etching → Microscopic determination

Microstructure images of specimens following for hardening zone and base material.

